

Product Datasheet

Human ZPI ELISA Kit (orb443000)

Description	Human Protein Z Dependent Protease Inhibitor (ZPI) Sandwich ELISA Kit
Reactivity	Human
Range	1.25-80ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	Q9UK55
Sensitivity	0.49ng/mL
Entrez	51156
Expiration Date	Please enquire.