

Product Datasheet

Mouse ZPI ELISA Kit (orb443001)

Description	Mouse Protein Z Dependent Protease Inhibitor (ZPI) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.937-60ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	Q8R121
Sensitivity	0.41ng/mL
Entrez	217847
Expiration Date	Please enquire.