

Product Datasheet

Rat STAT4 ELISA Kit (orb443018)

Description	Rat Signal Transducer And Activator Of Transcription 4 (STAT4) Sandwich ELISA Kit
Reactivity	Rat
Range	78.12-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q66HB2
Sensitivity	28.56pg/mL
Entrez	367264
Expiration Date	Please enquire.