

Product Datasheet

Equine INS ELISA Kit (orb443019)

Description	Equine Insulin (INS) competitive inhibition ELISA Kit
Reactivity	Equine
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P01310
Sensitivity	0.69ng/mL
Entrez	100060077
Expiration Date	Please enquire.