

Product Datasheet

Human TNFRSF8 ELISA Kit (orb443026)

Description	Human Tumor Necrosis Factor Receptor Superfamily, Member 8 (TNFRSF8) Sandwich ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P28908
Sensitivity	6.79pg/mL
Entrez	943
Expiration Date	Please enquire.