

Product Datasheet

Mouse DBN1 ELISA Kit (orb443036)

Description	Mouse Drebrin 1 (DBN1) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9QXS6
Sensitivity	6.3pg/mL
Entrez	56320
Expiration Date	Please enquire.