

Product Datasheet

Porcine DMT1 ELISA Kit (orb443040)

Description	Porcine Divalent Metal Transporter 1 (DMT1) Sandwich ELISA Kit
Reactivity	Porcine
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	B2ZDZ3
Sensitivity	0.062ng/mL
Entrez	733595
Expiration Date	Please enquire.