

Product Datasheet

Human BAI1 ELISA Kit (orb443058)

Description	Human Brain Specific Angiogenesis Inhibitor 1 (BAI1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O14514
Sensitivity	0.065ng/mL
Entrez	575
Expiration Date	Please enquire.