

Product Datasheet

Human IRF2BP2 ELISA Kit (orb443067)

Description	Human Interferon Regulatory Factor 2 Binding Protein 2 (IRF2BP2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q7Z5L9
Sensitivity	0.066ng/mL
Entrez	359948
Expiration Date	Please enquire.