

Product Datasheet

Human KCNJ10 ELISA Kit (orb443068)

Description	Human Potassium Inwardly Rectifying Channel Subfamily J, Member 10 (KCNJ10) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P78508
Sensitivity	0.64ng/mL
Entrez	3766
Expiration Date	Please enquire.