

Product Datasheet

Doublecortin antibody (orb443079)

Description	Goat polyclonal antibody to Doublecortin
Species/Host	Goat
Reactivity	Bovine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC
Target	Doublecortin / DCX (aa69-81)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: Paraffin embedded Human Brain (Cortex). Recommended concentration: 6µg/ml.. Notes: Immunofluorescence: Strong expression of the protein seen in the cytoplasm of HepG2 cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of Kelly cells. Recommended concentration: 10ug/m
Protein Sequence	DRSHCQSLRFHQN
Clonality	Polyclonal
MW	49.3
Entrez	1641
Dilution Range	P-ELISA: 1: 16000, IF/ICC: 10µg/ml, FACS: 10 ug/ml
Expiration Date	12 months from date of receipt.

