

Product Datasheet

Cadherin 5 antibody (orb443083)

Description	Goat polyclonal antibody to Cadherin 5
Species/Host	Goat
Reactivity	Bovine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Target	VE-cadherin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Notes: Immunofluorescence: Strong expression of the protein seen in the cytoplasm of HeLa cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of HeLa cells. Recommended concentration: 10ug/ml.
Protein Sequence	KQARAHGKSVPE
Clonality	Polyclonal
MW	87.5
Entrez	1003
Dilution Range	P-ELISA: 1: 128000, IF/ICC: 10µg/ml, FACS: 10µg/ml
Expiration Date	12 months from date of receipt.