



Product Datasheet

HYAL1 antibody (orb443084)


Description	Goat polyclonal antibody to HYAL1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC
Target	Hyaluronidase 1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: Paraffin embedded Human Liver. Recommended concentration: 5µg/ml.. Notes: Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A431 and HeLa cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of HepG2 cells. Recommended concentration:
Protein Sequence	ENTRTKESCQAIKEY
Clonality	Polyclonal
MW	48.4; 45.1; 28.0; 19.4
Entrez	3373
Dilution Range	ELISA: 1:128000, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of HepG2 cells ...



Immunofluorescence analysis of A431 cell...



Immunofluorescence analysis of HeLa cell...