



Product Datasheet

CD5 Antibody (orb443097)


Description	CD5 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, IHC, IHC-Fr, WB
Immunogen	E. coli-derived human CD5 recombinant protein (Position: R25-L495).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	67 kDa
Uniprot ID	P06127
Dilution Range	Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml Immunohistochemistry (Frozen Section), 0.5-1µg/ml Immunocytochemistry, 0.5-1µg/ml Immunofluorescence, 2µg/ml Flow Cytometry, 1-3µg/1x10 ⁶ cells Direct ELISA, 0.1-0.5µg/ml
Expiration Date	12 months from date of receipt.



Flow Cytometry analysis of EL-4 cells us...



WB analysis of CD5 using anti-CD5 antibo...



IF analysis of CD68 and CD5 using anti-C...