



Product Datasheet

Cytochrome P450 17A1/Cyp17a1 Antibody (orb443105)


Description	Cytochrome P450 17A1/Cyp17a1 Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	E. coli-derived mouse Cyp17a1 recombinant protein (Position: Q80-R363).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57 kDa
Uniprot ID	P27786
Dilution Range	Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml Immunofluorescence, 2µg/ml Direct ELISA, 0.1-0.5µg/ml
Expiration Date	12 months from date of receipt.



WB analysis of Cyp17a1 using anti-Cyp17a...



IF analysis of Cyp17a1 using anti-Cyp17a...



IF analysis of Cyp17a1 using anti-Cyp17a...