



## Product Datasheet

### IL36 alpha/IL36A Antibody (orb443114)

|                            |                                                                                                                                                                                                       |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | IL36 alpha/IL36A Antibody                                                                                                                                                                             |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                          |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                                                                                                                                                        |
| <b>Immunogen</b>           | E. coli-derived human IL36 alpha recombinant protein (Position: K6-D151).                                                                                                                             |
| <b>Form/Appearance</b>     | Lyophilized                                                                                                                                                                                           |
| <b>Concentration</b>       | Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.                                                                                                                             |
| <b>Storage</b>             | Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                                                                                 |
| <b>Application notes</b>   | Add 0.2ml of distilled water will yield a concentration of 500ug/ml.                                                                                                                                  |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                                            |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                            |
| <b>MW</b>                  | 25 kDa                                                                                                                                                                                                |
| <b>Uniprot ID</b>          | <a href="#">Q9UHA7</a>                                                                                                                                                                                |
| <b>Dilution Range</b>      | Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml ELISA (Cap), 1-5µg/ml                                                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                       |



WB analysis of IL36 alpha using anti-IL36...



IHC analysis of IL36 alpha using anti-IL36...