


Product Datasheet

Human AGT/Serpin A8/Angiotensinogen ELISA Kit (orb443235)

Description	Human AGT/Serpin A8/Angiotensinogen ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D34-A485
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (EDTA).
Uniprot ID	P01019
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:500, concentration in serum and plasma is 1-2 ug/ml.
Expiration Date	Please enquire.



Human
AGT
PicoKine
ELISA Kit
standard
cu...