


Product Datasheet

Human Cystatin SN/CST1 ELISA Kit (orb443246)

Description	Human Cystatin SN/CST1 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: W21-S141
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva and urine.
Uniprot ID	P01037
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is around 2 ng/ml.
Expiration Date	Please enquire.



Human
Cystatin
SN/CST1
PicoKine
ELISA Ki...