


Product Datasheet

Human Autotaxin ELISA Kit (orb443247)

Description	Human Autotaxin ELISA Kit
Reactivity	Human
Range	625 pg/ml - 40,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A36-I863
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin) , urine and human milk.
Uniprot ID	Q13822
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:20, concentration in serum is around 200ng/ml
Expiration Date	Please enquire.



Human
Autotaxin
PicoKine
ELISA Kit
stand...