


Product Datasheet

Human IL-38/IL1F10 ELISA Kit (orb443253)

Description	Human IL-38/IL1F10 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: M1-W152
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q8WWZ1
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 62.5 pg/ml.
Expiration Date	Please enquire.



Human IL-38/IL1F10 PicoKine ELISA Kit st...