


Product Datasheet

Human MCAM/CD146 ELISA Kit (orb443259)

Description	Human MCAM/CD146 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: V24-G559
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
Uniprot ID	P43121
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is around 1500 pg/ml.
Expiration Date	Please enquire.



Human
MCAM/CD146
PicoKine
ELISA Kit
stan...