


Product Datasheet

Rat TFF3/ITF ELISA Kit (orb443260)

Description	Rat TFF3/ITF ELISA Kit
Reactivity	Rat
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A23-F81
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q03191
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 100 pg/ml.
Expiration Date	Please enquire.



Rat TFF3
PicoKine
ELISA Kit
standard
cur...