

Product Datasheet

CD43 antibody (orb44350)

Description	Mouse Monoclonal to CD43.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, IP, WB
Immunogen	Human T lymphocytes.
Target	CD43
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1 µg/ml. Immunohistochemistry (paraffin sections): Positive tissue: spleen, thymus, tonsil; recommended dilution: 10 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P16150
Entrez	6693
Dilution Range	Flow cytometry: Recommended dilution: 1 µg/ml. Immunohistochemistry (paraffin sections): Positive tissue: spleen, thymus, tonsil; recommended dilution: 10 µg/ml.
Expiration Date	12 months from date of receipt.

