



Product Datasheet

HLA-A antibody (orb44401)

Description	Mouse monoclonal antibody to HLA-A
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Leukocytes of a patient suffering with LGL type leukemia.
Target	HLA Class I
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Dilution Range	Flow cytometry: Recommended dilution: 1 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human leukocytes stained u...



Flow cytometry surface staining pattern ...