

Product Datasheet

IFN gamma antibody (orb444035)

Description	Mouse monoclonal antibody to IFN gamma
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF
Concentration	See vial for concentration
Note	For research use only
Application notes	Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P01579
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.