


Product Datasheet

CD8 antibody (orb44413)

Description	Mouse monoclonal antibody to CD8
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	Leukocytes of a patient suffering from LGL Type Leukemia.
Target	CD8
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2 µg/ml. Immunoprecipitation: Excellent for immunoisolation of CD8 positive T cells.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Dilution Range	Flow cytometry: Recommended dilution: 2 µg/ml. Immunoprecipitation: Excellent for immunoisolation of CD8 positive T cells.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...