



Product Datasheet

HLA-E antibody (orb44422)


Description	Mouse Monoclonal to HLA-E.
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, IP
Immunogen	Bacterially expressed recombinant HLA-E refolded with beta2-microglobulin and peptide.
Target	HLA-E
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml; positive tissue: spleen. Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P13747
Entrez	3133
Dilution Range	Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml; positive tissue: spleen. Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of human monocytes (red-filled...



Flow cytometry surface staining pattern ...



Immunohistochemistry staining of human s...