

Product Datasheet

HLA-G antibody (orb44426)

Description

Mouse Monoclonal to HLA-G.

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-Fr, IHC-P, WB

Applications

Immunogen

Denatured bacterially expressed recombinant human HLA-G heavy chain.

Target

HLA-G

Preservatives

Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Concentration

1 mg/ml

Storage

Store at 2-8°C. Do not freeze.

Note

For research use only

Application notes

Immunohistochemistry (paraffin sections):
Recommended dilution: 1:60 to 1:100, incubation: 1 h at RT; positive tissue: human placenta - extravillous cytotrophoblast, heat retrieval in 0.01M citrate buffer (4x2 min. in microwave oven).

Isotype

Mouse IgG1

Clonality

Monoclonal

Purity

Purified by protein-A affinity chromatography.

Uniprot ID

P17693

Entrez


3135

Dilution Range


Immunohistochemistry (paraffin sections):
Recommended dilution: 1:60 to 1:100, incubation: 1 h at RT; positive tissue: human placenta - extravillous cytotrophoblast, heat retrieval in 0.01M citrate buffer (4x2 min. in microwave oven).Western blotting:
Recommended dilution 1-2 µg/ml.

Expiration Date


12 months from date of receipt.



Immunohistochemistry staining with anti-...



Immunohistochemistry staining with anti-...



Western blotting analysis of human HLA-G...