


Product Datasheet

HLA-G antibody (PE) (orb44438)

Description	Mouse Monoclonal to HLA-G.
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	Recombinant human HLA-G refolded with beta2-microglobulin and peptide.
Target	HLA-G
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.1 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml; positive control: JEG-3 human choriocarcinoma cell line.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P17693
Entrez	3135
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml; positive control: JEG-3 human choriocarcinoma cell line.
Expiration Date	12 months from date of receipt.



Separation of
LCL 721.221
HLA-G
transfec...