



Product Datasheet

WRB antibody (orb44448)


Description	WRB antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Tail-anchored protein insertion receptor WRB protein (38-99AA)
Target	WRB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	000258
Dilution Range	WB:1:2000-5000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse spleen ti...



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of human te...