

Product Datasheet

C CBL (phospho-Tyr674) antibody (orb4446)

Description	C CBL (phospho-Tyr674) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human C CBL around the phosphorylation site of Tyr674 (Al(p-Y)SL)
Target	phospho-CBL(Tyr674)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100 kDa
Uniprot ID	P22681
Expiration Date	12 months from date of receipt.