

Product Datasheet

Neurofilament medium protein antibody (orb44467)

Description

Mouse Monoclonal to Neurofilament medium

Reactivity

Mammal

Conjugation

Unconjugated

Tested

ICC, IHC-P, WB

Applications

Immunogen

Pellet of porcine brain cold stable proteins after depolymerization of microtubules.

Target

Neurofilament M

Preservatives

Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Concentration

1 mg/ml

Storage

Store at 2-8°C. Do not freeze.

Note

For research use only

Application notes

Immunocytochemistry: Positive tissue: Neuro2A murine cell line, Carnoys fixative 2 x 3 min, blocking 1% glycine + 0.2% gelatin 10 min. Western blotting: Recommended dilution: 1-2 µg/ml.

Isotype

Mouse IgG2a

Clonality

Monoclonal

Purity

Purified by protein-A affinity chromatography.

Uniprot ID

[P07197](#)

Entrez


[4741](#)

Dilution Range


Immunocytochemistry: Positive tissue: Neuro2A murine cell line, Carnoys fixative 2 x 3 min, blocking 1% glycine + 0.2% gelatin 10 min. Western blotting: Recommended dilution: 1-2 µg/ml.

Expiration Date

12 months from date of receipt.



Immunocytochemistry staining of neurofil...



Western blotting analysis of human neuro...