

Product Datasheet

PAG antibody (orb44470)

Description

Mouse monoclonal antibody to PAG

kDa 1 2 3

Reactivity

Bovine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, IP, WB

Applications
Immunogen

C-terminal peptide (last 15 amino acids) of human Csk binding protein coupled to KLH.

Immunoprecipitation
of human PAG/Cbp
fro...

Target

PAG1

Preservatives

Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Concentration

1 mg/ml

Storage


Store at 2-8°C. Do not freeze.

Note

For research use only

Application notes

Immunohistochemistry (paraffin sections): Positive tissue: appendix (germinal center of lymphatic follicle), heat-mediated antigen retrieval in citrate buffer pH 6.1. Immunoprecipitation: Positive control: RAJI human Burkitt lymphoma cell line.


Immunohistochemistry
staining (paraffin ...

Isotype

Mouse IgG2b

Clonality

Monoclonal

Purity

Purified by protein-A affinity chromatography.

Uniprot ID
Q3U1F9
Entrez
94212
Dilution Range

Immunohistochemistry (paraffin sections): Positive tissue: appendix (germinal center of lymphatic follicle), heat-mediated antigen retrieval in citrate buffer pH 6.1. Immunoprecipitation: Positive control: RAJI human Burkitt lymphoma cell line.

Expiration Date

12 months from date of receipt.