


Product Datasheet

STAT1 (phospho-Ser727) antibody (orb44474)

Description	Mouse Monoclonal to STAT1 (phospho-Ser727).
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	STAT1 peptide sequence 721-733 (Ser727 phosphorylated).
Target	STAT1 (pS727)
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P42224
Entrez	6772
Expiration Date	12 months from date of receipt.



Western blotting analysis (reducing cond...)