


Product Datasheet

Transferrin antibody (orb44475)

Description	Mouse Monoclonal to Transferrin.
Reactivity	Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Porcine transferrin
Target	Transferrin
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P09571
Entrez	396996
Expiration Date	12 months from date of receipt.



Western blotting analysis (reducing cond...