



Product Datasheet

SHIP-1 antibody (orb44487)

Description	Mouse Monoclonal to SHIP-1.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Peptide corresponding to a sequence within N-terminal domain of Human SHIP-1.
Target	SHIP1
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Intracellular staining; recommended dilution: 2-5 µg/ml; positive control: human blood leukocytes. Western blotting: Positive control: RAMOS human cell line, reducing conditions.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q92835
Entrez	3635
Dilution Range	Flow cytometry: Intracellular staining; recommended dilution: 2-5 µg/ml; positive control: human blood leukocytes. Western blotting: Positive control: RAMOS human cell line, reducing conditions.
Expiration Date	12 months from date of receipt.



Western blotting analysis (reducing cond...



Separation of MOLT-4 cells stained using...