

Product Datasheet

Syk antibody (orb44501)

Description

Mouse Monoclonal to Syk.

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

Applications

ICC, IHC-P, IP, WB

Immunogen

Recombinant fragment (aa 5-360) of human Syk.

Target

SYK

Preservatives

Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Concentration

1 mg/ml

Storage

Store at 2-8°C. Do not freeze.

Note

For research use only

Application notes

Western blotting: Recommended dilution: 1-2 µg/ml; positive control: RBL rat basophilic leukemia cell line, A-431 human epidermoid carcinoma cell line, RAMOS lymphoma cell line, U-937 human histiocytic lymphoma cell line, JURKAT human peripheral blood T cell leukemia cell line; negative control: HeLa human cervix carcinoma cell line; non-reducing conditions. Immunohistochemistry (paraffin sections): Recommended dilution: 5 µg/ml; positive tissue: tonsil B cells.

Isotype

Mouse IgG1

Clonality

Monoclonal

Purity

Purified by protein-A affinity chromatography.

Uniprot ID

[P43405](#)

Entrez



[6850](#)

Dilution Range


Western blotting: Recommended dilution: 1-2 µg/ml. Immunohistochemistry (paraffin sections): Recommended dilution: 5 µg/ml; positive tissue: tonsil B cells.

Expiration Date


12 months from date of receipt.



Immunocytochemistry analysis of Syk in h...



Immunocytochemistry analysis of Syk in h...



Western blotting analysis of human Syk u...