

Product Datasheet

HSV1 antibody (orb44505)

Description	Mouse Monoclonal to HSV1.
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	Exctract of HSV-1 infected VERO (green monkey kidney) cells.
Target	Glycoprotein C (HSV-1)
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
Purity	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Uniprot ID	Q8UYE2
Entrez	2703410
Expiration Date	12 months from date of receipt.