


Product Datasheet

CD162 antibody (orb44511)

Description	Mouse Monoclonal to CD162.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human thymocytes
Target	CD162
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1.5 - 2.5 µg/ml; positive control: human peripheral blood.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q14242
Entrez	6404
Dilution Range	Flow cytometry: Recommended dilution: 1.5 - 2.5 µg/ml; positive control: human peripheral blood.
Expiration Date	12 months from date of receipt.



Flow cytometry (surface staining) of hum...