



Product Datasheet

CD36 antibody (orb44516)

Description	Mouse monoclonal antibody to CD36
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Platelets
Target	CD36
Preservatives	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P16671
Entrez	948
Dilution Range	Flow cytometry: Recommended dilution: 2 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of murine CD36 positive throm...



Flow cytometry surface staining pattern ...