



Product Datasheet

TRIM Antibody (PE) (orb44527)


Description	Anti-TRIM PE
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	Recombinant intracellular domain (aa 29-186) of human TRIM.
Target	TRIM
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.1 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	Q6PIZ9
Entrez	50852
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.
Expiration Date	12 months from date of receipt.



Flow cytometry multicolor surface staini...



Comparison of anti-TRIM (TRIM-04) PE ant...



Separation of TRIM positive CD3 positive...