

Product Datasheet

LSP1 (phospho-Ser252) antibody (orb4453)

Description	LSP1 (phospho-Ser252) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human LSP1 around the phosphorylation site of Ser252 (QA(P-S)IE)
Target	phospho-LSP1 (Ser252)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37 kDa
Uniprot ID	P33241
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.