



Product Datasheet

HLA-A antibody (Biotin) (orb44577)

Description	Mouse monoclonal antibody against HLA-A conjugated to Biotin
Reactivity	Bovine, Feline, Human, Primate
Conjugation	Biotin
Tested Applications	ELISA, FC, ICC, IHC-Fr, IP
Immunogen	Membrane of human tonsil cells
Target	HLA Class I
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-12 µg/ml
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Dilution Range	Flow cytometry: Recommended dilution: 1-12 µg/ml
Expiration Date	12 months from date of receipt.



Separation of human leukocytes (red-fill...)



Flow cytometry surface staining pattern ...