



Product Datasheet

HLA-A antibody (FITC) (orb44578)

Description	Mouse monoclonal antibody against HLA-A conjugated to FITC
Reactivity	Bovine, Feline, Human, Primate
Conjugation	FITC
Tested Applications	FC
Immunogen	Membrane of human tonsil cells
Target	HLA Class I
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-3 µg/ml.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Dilution Range	Flow cytometry: Recommended dilution: 1-3 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human leukocytes (red-fill...)



Flow cytometry surface staining pattern ...