



Product Datasheet

CD7 antibody (orb44604)

Description	Mouse Monoclonal to CD7.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, IP, WB
Immunogen	Human acute myelogenous leukaemia cell line KG-1.
Target	CD7
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Immunohistochemistry (paraffin sections): Recommended dilution: 5 µg/ml; positive tissue: spleen. Flow cytometry: Recommended dilution: 2 µg/ml. Western blotting: Non-reducing conditions.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P09564
Entrez	924
Dilution Range	Immunohistochemistry (paraffin sections): Recommended dilution: 5 µg/ml; positive tissue: spleen. Flow cytometry: Recommended dilution: 2 µg/ml. Western blotting: Non-reducing conditions.
Expiration Date	12 months from date of receipt.



Separation of human CD7 positive lymphoc...



Flow cytometry surface staining pattern ...