


Product Datasheet

FoxP3 antibody (PE) (orb44650)

Description	Mouse Monoclonal to FoxP3.
Reactivity	Human, Mouse
Conjugation	PE
Tested Applications	FC
Immunogen	Full-length His-tagged recombinant murine FoxP3
Target	FoxP3
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.5 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	Q9BZS1
Entrez	50943
Dilution Range	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(intracellular
s...