

Product Datasheet

RNA polymerase II CTD repeat YSPTSPS (phospho-Ser5) antibody (Biotin) [Out of stock] (orb446620)

Description

Rabbit polyclonal antibody against

YSPTSPS (phospho-Ser5)
conjugated to...

Species/Host

Rabbit

Conjugation

Biotin

Tested

ELISA, IF, IHC-Fr, IHC-P, WB

Applications
Immunogen

KLH conjugated Synthesised
phosphopeptide derived from
human RNA polymerase II CTD
repeat YSPTSPS around the
phosphorylation site of phospho
S5.

Preservatives

0.01M TBS(pH7.4) with 1%
rAlbumin, 0.03% Proclin300 and
50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid
repeated freeze/thaw cycles. The
lyophilized antibody is stable at
room temperature for at least one
month and for greater than a year
when kept at -20°C. When
reconstituted in sterile pH 7.4
0.01M PBS or diluent of antibody
the antibody is stable for at least
two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

240kDa

Uniprot ID
P08775, P24928
Dilution Range

WB: 1:50-200, ELISA: 1:100-1000,
IHC-P: 1:50-200, IF: 1:50-200

Expiration Date

12 months from date of receipt.