

Product Datasheet

CD32b antibody (FITC) (orb44665)

Description	CD32b antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Low affinity immunoglobulin gamma Fc region receptor II-b protein (12-289AA)
Target	FCGR2B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P31994
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.