


Product Datasheet

CD90 antibody (FITC) (orb44670)

Description

Mouse Monoclonal antibody against CD90 conjugated to



Separation of
JURKAT cells
(red) from
SP...

Reactivity

Equine, Human, Porcine, Primate

Conjugation

FITC

Tested

FC

Applications

Immunogen

HEL erythroleukemia cells

Target

CD90

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.

Isotype

Mouse IgG1 kappa

Clonality

Monoclonal

Purity

Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

P04216

Entrez

7070

Dilution Range

Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.

Expiration Date

12 months from date of receipt.