

Product Datasheet

CD324 antibody (orb44672)

Description	Mouse Monoclonal to CD324.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-Fr, IP, WB
Immunogen	T-47D cells
Target	CD324
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Immunohistochemistry (frozen sections): Recommended dilution: 4-8 µg/ml; positive tissue: tonsil. Flow cytometry: Recommended dilution: 5-10 µg/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P12830
Entrez	999
Dilution Range	Immunohistochemistry (frozen sections): Recommended dilution: 4-8 µg/ml; positive tissue: tonsil. Flow cytometry: Recommended dilution: 5-10 µg/ml
Expiration Date	12 months from date of receipt.